



FOR IMMEDIATE HEMORRHAGE CONTROL



Clamp iT

Clamp and seal the wound
with minimal pain



Contain iT

Contained hematoma places
direct pressure on the
injured vessels



Control iT

Control your patient's
emergency in seconds


iTraumacare.com

 iTClamp[®]

To order the iTClamp^(R),
email us at info@itraumacare.com,
or visit iTraumaCare.com.

iTClamp^(R) SPECIFICATIONS

PART NUMBER	9100
SIZE	7.09 cm x 7.21 cm x 4.09 cm
WEIGHT	1.3 oz / 37.1g
MATERIAL	Medical Grade Polycarbonate, 304 Stainless Steel
OTHER	Sterile, Latex Free



HEMORRHAGE CONTROL SYSTEM

 **iTClamp^(R)**
by INNOVATIVE
TRAUMA CARE^(R)